

## FUME HOOD

### ILB - 159 FUME HOOD

Fume Hood is fabricated of Stainless Steel/Mild Steel. It is designed so as to throw out all toxic/harmful vapors & fumes thus protecting costly instruments and the persons working in the lab. Working table top is acid/alkali resistant and a small wash basin with connections provided for inlet and outlet of water. The front door moves vertically up and down with concealed counter balanced weights. The unit is fitted with fluorescent light and a gas cock for gas/air supply. The front facing panel is fitted with 15/5 Amp socket with switches for exhaust system and fluorescent light. To work on 220 volts A.C.50 Hz supply.

#### Optional :

- \* High CFM exhaust System (for duct exceeding 10 feet in length)
- \* S.S Lining on work bench.
- \* PVC Fume Duct (per ft.)
- \* Digital Controller
- \* Marble on Working Table
- \* 3 Way Nozzels

#### Also Available in :

- \* Complete Stainless Steel (GMP Model)
- \* Complete inner outer PP Model

	A	B	C	D
Working Chamber	2'x2'x2'	3'x2'x2'	4'x2'x2'	6'x2'x2'



**NOTE : CUSTOMISED SIZES CAN BE MADE ACCORDING TO THE REQUIREMENT**

Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

OFFICE : # 17-18 SHREE RADHEY COMPLEX, NEAR B.D. HIGH SCHOOL, SCIENCE MARKET, AMBALA CANTT.

WORKS : # 57 VIKAS PURI, ADJACENT INDUSTRIAL AREA, SAI BABA MANDIR ROAD, AMBALA CANTT.

0171 - 4000114, +91 8295332770, +91 9896575277, +91 9729672661

Email : indolab.india@gmail.com - taman.indolab@gmail.com

www.indolabinstrument.com - www.indolab.net