

HUMIDITY CHAMBER

ILB-149 HUMIDITY CHAMBER

Standard Double walled construction Inner chamber made out of highly polished stainless steel sheet of 304/202 grade and outer made of mild sheet duly powder coated, temperature range 5°C above ambient to 60°C and humidity from 40% to 90%. Humidity & Temperature is controlled by PID type digital controller.

OPTIONAL :

- * Low water level cut off device.
- * Temperature and humidity overshoot alarm.
- * Microprocessor Based Digital Temperature Controller.
- * Automatic electronic controlled With Digital Display for Temperature & Humidity combined.
- * Automatic Voltage Stabilizer

Capacity in Litres

95 Litres

147 Litres

225 Litres

336 Litres

420 Litres



**NOTE : CUSTOMISED SIZES CAN BE MADE ACCORDING TO THE REQUIREMENT
COMPLETE STAINLESS STEEL (GMP) MODELS ARE ALSO AVAILABLE**

Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

OFFICE : # 17-18 SHREE RADHEY COMPLEX, NEAR B.D. HIGH SCHOOL, SCIENCE MARKET, AMBALA CANTT.

WORKS : # 57 VIKAS PURI, ADJACENT INDUSTRIAL AREA, SAI BABA MANDIR ROAD, AMBALA CANTT.

0171 - 4000114, +91 8295332770, +91 9896575277, +91 9729672661

Email : indolab.india@gmail.com - taman.indolab@gmail.com

www.indolabinstruments.com - www.indolab.net