

ILB-107 BACTERIOLOGICAL INCUBATOR

Bacteriological incubator are environment friendly for bacteriological test and storage of sensitive cultures, vaccines, culture of bacterial micro organisms and plant life, serum incubation studies, Immunological work, chemicals etc.

FEATURES :

- * Bacteriological incubator is mainly used is drying of slides, tissue culture, incubation of anti body, paraffin embedding incubation of hydroxy steroids, crystalization and microbiological determination.
- * Bacteriological incubator is double walled construction. Inner chamber made of stainless steel.
- * Outer chamber is made of mild/stainless steel sheet, finished with bacteria resistant powder coating.
- * Special grade glass wool is used for insulation and to avoid heat losses.
- * Shelves are adjustable to set convenient height
- * Shelves are also made of stainless steel SS304/202.
- * The door has antifog double glass window for inspection of samples without opening the door.
- * Heater is made of high grade imported wire which are put inside the porcelain beads and placed at the bottom & both sides for uniform temperature all over the space.
- * Ventilation : The bacteriological incubator is fitted with air ventilators on both side to remove gasses and fumes.
- * Power requirement : 220V-50Hz, AC supply.



OPTIONAL :

- * Air Circulating Fan.
- * Digital PID Temperature Controller.
- * Microprocessor based Digital Temperature Controller with Timer.
- * Audio-Visual Alarm.
- * Timer 0-999 hours.

Chamber Size (H x W x D) in mm	Chamber Size (H x W x D) in inch	Capacity in Litres
300x300x300	12"x12"x12"	28 Litres
350x350x350	14"x14"x14"	45 Litres
400x400x400	16"x16"x16"	65 Litres
450x450x450	18"x18"x18"	95 Litres
600x450x450	24"x18"x18"	125 Litres
600x600x600	24"x24"x24"	224 Litres
900x600x450	36"x24"x18"	252 Litres
900x600x600	36"x24"x24"	336 Litres
900x900x900	36"x36"x36"	720 Litres



NOTE : CUSTOMISED SIZES CAN BE MADE ACCORDING TO THE REQUIREMENT

Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

OFFICE : # 18 SHREE RADHEY COMPLEX, NEAR B.D. HIGH SCHOOL, SCIENCE MARKET, AMBALA CANTT.

WORKS : # 57 VIKAS PURI, ADJACENT INDUSTRIAL AREA, SAI BABA MANDIR ROAD, AMBALA CANTT.

0171 - 4000114, +91 8295332770, +91 9896575277, +91 9729672661

Email : indolab.india@gmail.com - taman.indolab@gmail.com

www.indolabinstrument.com - www.indolab.net