

LAMINAR AIR FLOW CABINET

ILB-157 LAMINAR AIR FLOW (HORIZONTAL)

Outer body made of Mild Steel duly powder coated/Stainless Steel. Inner chamber made of Stainless Steel. Working table is made of Stain Finish S.S. Table top. Complete with HEPA filter & prefilters motor blower, static pressure manometer, cock for gas, air or vacuum line. Front & Side doors are made of plexiglass transparent thick sheet duly framed. Built in U.V. Tube Work table is illuminated by flourescent light. Two legs table top and two table supports are detachable & can be taken out & fixed within minutes, front door can be kept full open/half open, fully closed as desired, The equipment is fitted with safety attachments of international standards.

Optional :

- * Magnehelic Gauge for Pressure
- * Digital Controller

| | A | B | C | D |
|----------------------|-----------|-----------|----------|----------|
| Working Chamber | 2'x2'x2' | 3'x2'x2' | 4'x2'x2' | 6'x2'x2' |
| Size of HEPA Filter | 2'x2'x6" | 3'x2'x6" | 4'x2'x6" | 3'x2'x6" |
| No. of HEPA Filter | 1 | 1 | 1 | 2 |
| No. of Pre Filter | 1 | 1 | 1 | 2 |
| Illumination | 1 x 20W | 1 x 20W | 1 x 20W | 2 x 20W |
| U.V.Germicidal Light | 1 x 1½ ft | 1 x 1½ ft | 1 x 3 ft | 2 x 3 ft |

Also Available in :

- * Complete Stainless Steel (GMP Model)



Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

ILB - 157 A LAMINAR AIR FLOW CABINET VERTICAL

All exterior surfaces made of mild steel duly powder coated/Stainless Steel. Inner chamber made of Stainless Steel. Working table is made of Stain finish S.S. Table top. Complete with HEPA filter & prefilters motor blower, static pressure manometer, cock for gas, air or vacuum line. Front doors is Made of plexiglass transparent thick sheet duly framed. Built in U.V. tube work table is Illuminated by flourescent light. Whole LAF cabinet is fitted on detachable table of M.S. powder coated/enamel paint. Front door can be kept full open/half open, fully closed as desired. But with laminar flows down vertically from top towards the work table. Supplied complete with static finish stain finish s.s. table top. Transparent front door are made of 5mm plexi glass, static pressure manometer, built in U.V. germicidal light & flourescent light and Cock for gas, air or vacuum line. The equipment is fitted with safety attachments of international standards.

| | A | B | C | D | E |
|----------------------|-----------|-----------|----------|----------|----------|
| Working Chamber | 2'x2'x2' | 3'x2'x2' | 4'x2'x2' | 6'x2'x2' | 8'x2'x2' |
| Size of HEPA Filter | 2'x2'x6" | 3'x2'x6" | 4'x2'x6" | 3'x2'x6" | 4'x2'x6" |
| No. of HEPA Filter | 1 | 1 | 1 | 2 | 2 |
| No. of Pre Filter | 1 | 1 | 1 | 2 | 2 |
| Illumination | 1 x 20W | 1 x 20W | 1 x 20W | 2 x 20W | 2 x 40W |
| U.V.Germicidal Light | 1 x 1½ ft | 1 x 1½ ft | 1 x 3 ft | 2 x 3 ft | 2 x 3 ft |

Also Available in :

* Complete Stainless Steel (GMP Model)



Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

OFFICE : # 18 SHREE RADHEY COMPLEX, NEAR B.D. HIGH SCHOOL, SCIENCE MARKET, AMBALA CANTT.

WORKS : # 57 VIKAS PURI, ADJACENT INDUSTRIAL AREA, SAI BABA MANDIR ROAD, AMBALA CANTT.

0171 - 4000114, +91 8295332770, +91 9896575277, +91 9729672661

Email : indolab.india@gmail.com - taman.indolab@gmail.com

www.indolabinstruments.com - www.indolab.net