

HOT AIR OVEN

ILB-101 HOT AIR OVEN

DESCRIPTION :

Double walled, inner chamber made of non-rustable SS304/202, outer body made of MS duly powder coated. 3.0" gap between the walls is filled with glass/mineral wool. Three heating elements made from the best quality nickel chromium wire are used on refractory supports. The elements work at black heat and have a very long life. Covered heaters are placed under ribs at the bottom and sides. This arrangement ensures equal uniform distribution of heat. Temperature is controlled by Digital Temperature Indicator cum controller upto 250° C. Air circulating fan for uniformity of temperature. Fitted with air ventilators on the sides of the oven. Perforated SS shelves, 3 core wire and plug for use on 220 Volts A/C.

Optional :

- * Air Circulating Fan
- * Microprocessor based PID Digital Temperature, Indicator cum Controller with Timer
- * Complete stainless steel GMP models are also available
- * Thermal Printer



Available in sizes:

Chamber Size in mm	Chamber Size in inch	Capacity (approx.) Litres
300 x 300 x 300	12" x 12" x 12"	28
350 x 350 x 350	14" x 14" x 14"	45
400 x 400 x 400	16" x 16" x 16"	65
450 x 450 x 450	18" x 18" x 18"	95
600 x 450 x 450	24" x 18" x 18"	125
600 x 600 x 600	24" x 24" x 24"	225
900 x 600 x 450	36" x 24" x 18"	252
900 x 600 x 600	36" x 24" x 24"	336
900 x 900 x 900	36" x 36" x 36"	720
1200 x 900 x 900	48" x 36" x 36"	1020



NOTE : CUSTOMISED SIZES CAN BE MADE ACCORDING TO THE REQUIREMENT.

Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

OFFICE : # 18 SHREE RADHEY COMPLEX, NEAR B.D. HIGH SCHOOL, SCIENCE MARKET, AMBALA CANTT.

WORKS : # 57 VIKAS PURI, ADJACENT INDUSTRIAL AREA, SAI BABA MANDIR ROAD, AMBALA CANTT.

0171 - 4000114, +91 8295332770, +91 9896575277, +91 9729672661

Email : indolab.india@gmail.com - taman.indolab@gmail.com

www.indolabinstrument.com - www.indolab.net