

### ILB-150 ENVIRONMENTAL CHAMBER / STABILITY CHAMBER

Standard Double walled construction Inner chamber made out of highly polished stainless steel sheet of 304/202 grade and outer made of mild sheet duly powder coated, temperature range 10°C to 60°C and humidity from 35% to 85%, 5% or less Rh controlled by PID type digital controller.

#### FEATURES :

- \* Digital LCD display.
- \* Microscope based PID temperature controller cum indicator.
- \* Auto dehumidifier for better accuracy.
- \* Temperature below ambient is controlled by a sealed compressor. stainless steel water reservoirs heated by immersion heaters provide humidity complete with all accessories.

#### OPTIONAL :

- \* Low water level cut off device.
- \* Temperature and humidity overshoot alarm.
- \* Microprocessor Based Digital Temperature Controller.
- \* Automatic electronic controlled With Digital Display for Temperature & Humidity combined.
- \* Automatic Voltage Stabilizer



Capacity in Cubic Feet	Capacity in Litres
4	112 Litres
6	171 Litres
10	280 Litres
12	336 Litres
15	420 Litres

**NOTE : CUSTOMISED SIZES CAN BE MADE ACCORDING TO THE REQUIREMENT**

#### Manufacturers & Suppliers :

- Scientific Equipment • Pharmaceutical Equipment • Laboratory Instruments • Hospital Equipment • Surgical Equipment
- Agricultural Equipment • Dairy Equipment • Clean Air Equipment • School Laboratories

OFFICE : # 18 SHREE RADHEY COMPLEX, NEAR B.D. HIGH SCHOOL, SCIENCE MARKET, AMBALA CANTT.

WORKS : # 57 VIKAS PURI, ADJACENT INDUSTRIAL AREA, SAI BABA MANDIR ROAD, AMBALA CANTT.

0171 - 4000114, +91 8295332770, +91 9896575277, +91 9729672661

Email : [indolab.india@gmail.com](mailto:indolab.india@gmail.com) - [taman.indolab@gmail.com](mailto:taman.indolab@gmail.com)

[www.indolabinstrument.com](http://www.indolabinstrument.com) - [www.indolab.net](http://www.indolab.net)